


Dental Lab Teeth Feldspathic Porcelain Veneers Custom Shaped

Our Product Introduction

for more products please visit us on adschinadentalab.com

Basic Information




Product Specification

- Turnaround Time: 3 Days In Lab
- Warranty: 5 Years
- Shipping Way: UPS (preferred), Fedex, DHL, Linex, Armex
- Payment: T/T, HSBC, PayPal, MoneyGram, Western Union
- Highlight: Teeth feldspathic porcelain, Custom Shaped feldspathic porcelain, Custom Shaped feldspathic veneers



More Images



Product Description

Description

Introduction:

As a manufacturer located in China, we are proud to introduce our latest product – Feldspathic Porcelain Veneers. We specialize in wholesale products to merchants located in countries outside of China. Our company is known for its high-quality products, and Feldspathic Porcelain Veneers are no exception. In this product introduction, we will detail what makes our veneers unique and why they should catch the attention of merchants worldwide.

What are Feldspathic Porcelain Veneers?

Feldspathic Porcelain Veneers are a type of dental veneer. They are made from layers of thin, custom-shaped porcelain that mimic the natural color and translucency of teeth. These veneers are bonded to the front surface of your teeth, creating a natural-looking and long-lasting smile. Unlike composite veneers, Feldspathic Porcelain Veneers are handcrafted by experienced technicians, providing a level of customization and quality that is unmatched in the industry.

What Makes Our Feldspathic Porcelain Veneers Unique?

1. Handcrafted Excellence

Our Feldspathic Porcelain Veneers are created by highly skilled technicians who have years of experience in the industry. The process of handcrafting each veneer ensures that every detail is perfect, from the shape, color, and translucency to the texture.

and finish. This level of customization is something that is not available with other types of veneers.

2. High-Quality Materials

The materials used to create our Feldspathic Porcelain Veneers are of the highest quality. We use the finest porcelain available, which is known for its strength, durability, and natural appearance. Our materials have been tested and approved for dental use, ensuring that our products meet the highest standards of quality and safety.

3. Outstanding Aesthetics

The aesthetics of our Feldspathic Porcelain Veneers are unmatched. We carefully select and match the shade of the veneers to the natural color of the patient's teeth. The layers of porcelain are carefully shaped and stacked to mimic the natural thickness and translucency of the teeth, resulting in a natural-looking and beautiful smile.

4. Easy to Apply

Our Feldspathic Porcelain Veneers are easy to apply and require minimal tooth preparation. We aim to provide a pain-free and comfortable experience for patients while still achieving excellent results. Our veneers are bonded to the teeth using a specialized adhesive, ensuring a secure and long-lasting bond.

5. Long-Lasting Results

Our Feldspathic Porcelain Veneers provide long-lasting results. With proper care and maintenance, they can last up to 15 years or more. Our veneers are also resistant to staining and discoloration, ensuring that patients can maintain their beautiful smile for years to come.

Conclusion:

At our company, we pride ourselves on providing high-quality products at an affordable price. Our Feldspathic Porcelain Veneers are no exception. We believe that our handcrafted excellence, high-quality materials, outstanding aesthetics, easy application, and long-lasting results are what sets us apart from other manufacturers. We invite merchants worldwide to try our product and experience the quality and craftsmanship that our company is known for.

Production Introduction:

Feldspathic veneers have been made popular by patients for their non-invasive preparation requirements, and their ability to allow the existing tooth shades to shine through. This hand-crafted veneer is created with artistry by use of a brush and glass-based powder mixed with liquid. To fabricate feldspathic veneers, the technician hand paints the veneer layer by layer so contour and texture can be controlled and crafted. The technicians at China Digital Lab have the experience necessary to craft highly esthetic veneers that will fit every patient's cosmetic needs and desires. Our veneers are excellent cosmetic solutions in cases of chipped, discolored, or gapped teeth.

Product Parameter

Customized

Product Feature and Application

Feldspathic porcelain veneers are a highly esthetic restoration that is ideal for anterior restorations where shape, texture, and natural-beauty are priorities.

ADVANTAGES/BENEFITS

- Minimal Prep (0.2mm-0.3mm)
- Highly esthetic
- Non-invasive
- Translucent

This veneering fabrication style may be recommended for patients when:

- The drilled tooth still has most of its enamel
- The patient does not grind their teeth
- When shade matching is difficult
- With extremely demanding patients


Deliver,Shipping and Serving

We accept STL files.

We accept normal impressions, you can send us by Fedex,UPS,DHL or Linex

We finish the jobs in 4-5 days once we receive your cases

We accept TT.Paypal



Why choose us?

Our Dental Veneers are easy to apply and can transform the look and feel of your smile in no time.

We look forward to cooperating with more customers or business dealings, and work together for the development of the industry.

We take pride in our strict quality control processes to ensure the excellence of our Dental Veneers products.

We combine the actual development of the enterprise, with a pragmatic attitude and down-to-earth style, we constantly reform, reorganize and innovate the property rights system, management system and operation system of the enterprise.

We pride ourselves on providing exceptional customer service and after-sales support for our Dental Veneers products.

Every initiative of our company introduces the concept of environmental protection, health and safety to create more opportunities for our sustainable development.

Our Dental Veneers products are ideal for those looking for a cost-effective and non-invasive dental solution.

We have strong strength and rich experience in the Feldspathic Porcelain Veneers industry.

We offer fast and reliable shipping options for our Dental Veneers products, worldwide.

The service tenet of creating value for customers is the core of our company's continuous healthy and stable development.

Hot Tags: feldspathic porcelain veneers, China feldspathic porcelain veneers manufacturers, suppliers, factory, dental expanders for adults, expansion appliance in orthodontics, expansion plate dental, Feldspathic Porcelain Veneers, maxillary palatal expander, expansion appliance dental

ADS DENTAL LAB 卓艺美学中心 ADS Dental Lab

+8613430498390

info@chinaadsdentalab.com

adschinadentalab.com

Building 4, Hengguangyao Industrial Park, Yonghe Road, Fuyong, Shenzhen, China 518103