

Removable Telescopic Partial Denture China Dental Lab For Missing Teeth

Our Product Introduction


for more products please visit us on adschinadentallab.com

Basic Information





Product Specification

- Turnaround Time: 3 Days In Lab
- Warranty: 5 Years
- Shipping Way: UPS (preferred),Fedex,DHL Linex,Armex
- Payment: T/T,HSBC,PayPal,MoneyGram,Western Union
- Highlight: Removable telescopic partial denture,
Missing Teeth telescopic partial denture,
Lab telescopic removable partial denture



More Images



Product Description

Description

Telescopic Partial Dentures - A Revolutionary Solution for Missing Teeth

As a leading manufacturer of dental products located in China, we take pride in offering high-quality and innovative solutions to meet the needs of our customers. We understand the importance of having healthy teeth and our products are designed to ensure that our customers achieve optimal oral health.

We are excited to introduce our latest product, Telescopic Partial Dentures, which is a game-changer in the world of dentistry. This product is designed to provide a comfortable, natural-looking, and reliable solution for people who have missing teeth.

What are Telescopic Partial Dentures?

Telescopic Partial Dentures are custom-made dental appliances that consist of a metal framework which fits around the remaining natural teeth and a removable set of artificial teeth that snaps onto the metal frame. The metal frame has grooves and channels that fit into each other, creating a secure and stable foundation for the artificial teeth.

One of the unique features of Telescopic Partial Dentures is their telescopic mechanism, which allows the metal frame and the artificial teeth to slide in and out of each other, ensuring a perfect fit and maximum comfort.

Why Choose Telescopic Partial Dentures?

There are many reasons why Telescopic Partial Dentures are the preferred solution for people who have missing teeth.

1. Comfortable and Natural-Looking

Telescopic Partial Dentures are designed to fit snugly around the remaining natural teeth, creating a seamless and natural look. The telescopic mechanism ensures that the dentures are comfortable to wear, with no slipping or shifting.

2. Reliable and Durable

The metal framework of Telescopic Partial Dentures is made of high-quality materials, ensuring durability and longevity. The teeth are also made of high-quality materials, ensuring that they are resistant to wear and tear.

3. Customizable and Versatile

Telescopic Partial Dentures can be customized to fit the unique shape and size of each patient's mouth. They can also be used for patients who have missing teeth in different areas of the mouth, making them a versatile and flexible solution.

4. Easy to Maintain

Telescopic Partial Dentures are easy to maintain, requiring only regular cleaning and proper storage. They are also easy to remove and insert, allowing patients to take them out when needed.

5. Affordable

Compared to other dental solutions, Telescopic Partial Dentures are an affordable option for patients who want to restore their missing teeth.

Why Choose Us?

Our Telescopic Partial Dentures are made using state-of-the-art technology and equipment, ensuring that they are of the highest quality. We have a team of highly skilled and experienced technicians who are dedicated to providing the best possible service to our customers.

We adhere to strict quality control measures at every stage of the manufacturing process, ensuring that our products meet or exceed international standards. We are committed to providing our customers with the best possible products and services at affordable prices.

Conclusion

If you are looking for a reliable, comfortable, and natural-looking solution for missing teeth, then Telescopic Partial Dentures are the perfect choice. They offer a customizable, versatile, and affordable solution that is easy to maintain and provides long-lasting results.

Contact us today to learn more about our Telescopic Partial Dentures and how we can help you provide your patients with the best possible dental solutions.

Production Introduction:

The telescopic crown consists of an inner crown (also known as a primary or inner telescopic) and an outer crown (also known as a secondary or outer telescopic). It is used to attach removable dentures. The remaining teeth are used for anchoring. In conventional technique, telescopic crowns - also known as double crowns - were usually made of gold alloys or with the help of electroplating technology. With the advent of digital technology, new possibilities have opened up in the field of telescopic crown manufacture, such as zirconium, non-precious metal alloys, titanium or PEEK plastics.

Product Parameter

Customized

Product Feature and Application

Telescopic Partial Dentures have been successfully used in partially edentulous patients. This type of retainer provides guidance, support, and protection from dislodging movement, and it transfers bite forces along the long axis of the abutment teeth. Furthermore, the denture can easily be retrieved when an abutment needs to be extracted. Removing the prosthesis also makes cleaning easier. In contrast to other dental restorations, the prosthesis can be used for a long time. In the upper jaw, the telescopic prosthesis can in some cases avoid a palate plate, which make patients feel more pleasant. Since this type of restoration is completely coated with resin, it comes aesthetically very close to the natural dentition.

Deliver,Shipping and Serving


We accept STL files.

We accept normal impressions, you can send us by Fedex,UPS,DHL or Linex

We finish the jobs in 4-5 days once we receive your cases

We accept TT.Paypal





Why choose us?

We are dedicated to providing our clients with personalized attention and solutions.

More and more companies, including us, have clearly realized that Telescopic Partial Dentures containing green technology are the only asset for corporate survival and market competition.

Our dentures are designed to be easy to use and care for.

We always adhere to the high degree of unity of economic responsibility and social responsibility, and we do not forget to fulfill social responsibility while developing ourselves.

Our company is based in China and specializes in producing high-quality dental dentures.

We always achieve market competitiveness with innovation, pay attention to cultural construction, create a harmonious and enterprising atmosphere, and build a platform for the common development of individuals and the company.

Our focus is on delivering dental solutions that are both effective and aesthetically pleasing.

In order to improve the quality and safety of products, the company starts from the source and strictly monitors and manages production.

We are committed to using only the best materials and equipment for our dentures.

Meeting the market demand and producing marketable Telescopic Partial Dentures is one of the important issues we are facing at present.

Hot Tags: telescopic partial dentures, China telescopic partial dentures manufacturers, suppliers, factory, mouth appliance to widen jaw, hybrid zirconia bridge, Removable Acrylic Partial Denture, jaw expansion orthodontics, removable hawley wire retainers, chipped pfm crown

ADS DENTAL LAB 卓艺美学中心 **ADS Dental Lab**

+8613430498390

info@chinaadsdentallab.com

adschinadentallab.com

Building 4, Hengguangyao Industrial Park, Yonghe Road, Fuyong, Shenzhen, China 518103