

ODM Modified Hawley Removable Retainer Post Orthodontic Appliance

Our Product Introduction

for more products please visit us on adschinadentallab.com

Basic Information




Product Specification

- Turnaround Time: 3 Days In Lab
- Warranty: 5 Years
- Shipping Way: UPS (preferred), Fedex And DHL Linex. Armex Are Also Available
- Payment: T/T, HSBC, PayPal, MoneyGram, Western Union
- Highlight: **Modified Hawley Removable Retainer, Post Orthodontic Hawley Removable Retainer, ODM modified hawley retainer**



More Images



Product Description

Description

Introducing the Hawley Removable Retainer

Are you in search of a quality dental retainer that not only offers exceptional performance but also deliver remarkable comfort? Look no further than the Hawley Removable Retainer. This cutting-edge product is made using the latest technology to guarantee the best dental health outcomes.

The Hawley Removable Retainer is a removable orthodontic appliance that is designed to keep the teeth in their final position. It is typically used after orthodontic treatment to maintain the achieved alignment of teeth. This easy-to-use retainer is ideal for both post-orthodontic use and for individuals who want to maintain their dental alignment.

Here are some of the most important aspects of the Hawley Removable Retainer and how they can benefit your customers:

1. High-Quality Materials

The Hawley Removable Retainer is made using high-quality materials that are durable and easy to maintain. The retainer's design is created using acrylic material, while the wires that hold the appliance in place are made using stainless steel. These materials are not only strong but also highly resistant to wear and tear, guaranteeing longevity and maximum performance throughout their lifespan. This unique combination of quality materials allows the retainer to withstand daily use without significant wear or damage.

2. Removable Design

The Hawley Removable Retainer is designed to be removable quickly and easily, which makes it highly convenient and user-friendly. The wearer can easily take them out before meals, cleaning, or before any other oral activity, which makes it the ideal product for people with a busy lifestyle. The removable design also ensures that the retainer does not interfere with the user's daily activities, such as eating or speaking. Users can easily remove them, brush, floss and clean the retainer and put it back quickly and comfortably.

3. Comfortable Fit

The Hawley Removable Retainer is designed to fit flawlessly and comfortably within the mouth. It is a custom-made product that is fabricated to fit the user's specific teeth alignment perfectly, which ensures a comfortable and snug fit. This customized fit also ensures maximum effectiveness in maintaining teeth alignment and promoting good oral hygiene.

4. Aesthetic Appeal

The Hawley Removable Retainer comes in a variety of colors and designs that can blend seamlessly with the user's natural teeth alignment to match their unique style. The retainer can be designed to complement the user's dental aesthetics, offering them an effective and aesthetically-pleasing solution.

5. Easy Replacement

The Hawley Removable Retainer is designed to be easy to replace in case of damage or loss. We understand that accidents happen, which is why our retainer is designed for easy replacement. With our superior customer service, our clients can rest assured that we will provide outstanding support to promptly replace any worn or damaged retaining appliances.

In conclusion, the Hawley Removable Retainer is a high-quality orthodontic appliance that offers exceptional value to clients who want to maintain their dental alignment or after orthodontic treatment. It is made using materials that guarantee longevity, easy maintenance, and optimal performance. The removable design ensures it is convenient and comfortable to wear, while the customized fit, aesthetic appeal, and easy replacement make it the ideal solution for recovering patient's continued dental health. Whether you are a dental practitioner, dental distributor, or independent retailer, you can trust in Hawley Retainers and their excellent products' dedication to outstanding customer satisfaction. Contact us today for more information on this remarkable product.

Production Introduction:

A Hawley retainer is an oral appliance made of wire and hard plastic or acrylic. The wire extends across the teeth to keep them aligned after braces are removed. The acrylic rests on the roof of the mouth.

Product Parameter

Customized

Product Feature and Application

Hawley retainer is a removable orthodontic appliance and caring for this removable device is paramount for successful retention of your dream smile.

Caring for Your Hawley Retainer

- Cleaning daily to remove bacteria and food debris.

- Keeping away from heat or hot water.
- Storing in its case when not in the mouth.

How to Clean Hawley Retainer?

The Hawley retainer is worn during the night. The retainer should be washed with warm water and scrubbed lightly with a toothbrush. This is to be done before placing it in storage, and the storage for the container should also be rinsed regularly.

Deliver,Shipping and Serving

We accept STL files.

We accept normal impressions, you can send us by Fedex,UPS,DHL or Linex

We finish the jobs in 4-5 days once we receive your cases

We accept TT.Paypal





Why choose us?

Our Orthodontic Retainers products are backed by a comprehensive warranty for added peace of mind.

We value efficient coordination in all aspects, including corporate culture management, brand promotion, marketization and international competitiveness.

Our team of experts produces Orthodontic Retainers using state-of-the-art technology and advanced equipment.

The company continues to be user-experience-oriented and grandly launches self-developed products.

Our Orthodontic Retainers products are designed to be easy to use and operate for maximum convenience.

We have undertaken national, provincial and ministerial-level scientific research projects for many years, with fruitful technological achievements and outstanding strength.

We prioritize customer satisfaction above all else when it comes to our Orthodontic Retainers products.

Our enterprise spirit is the advanced group consciousness formed for a long time by the enterprise in production and operation activities to seek its own survival and development, and recognized by the masses of workers.

We are committed to delivering our Orthodontic Retainers products on time and in perfect condition.

We are well-established and professionally managed company. Being a leading manufacturer and supplier of Hawley Removable Retainer, we cater to various national and international clients for several years. With an effective and modern production capacity, we can fulfilled the varied requirements of clients.

Hot Tags: hawley removable retainer, China hawley removable retainer manufacturers, suppliers, factory, hawley retainer with bite plate, expansion in orthodontics, bottom teeth expander, dental arch expansion, cleaning palate expander, orthodontic expander for adults

ADS DENTAL LAB 卓艺美学中心 **ADS Dental Lab**

+8613430498390

info@chinaadsdentallab.com

adschinadentallab.com

Building 4, Hengguangyao Industrial Park, Yonghe Road, Fuyong, Shenzhen, China 518103