


Metal Cylindrical Precision Attachment Denture Telescoper Crown For Dental Prosthesis

Basic Information



Product Specification

- Shape: Cylindrical
- Usage: Dental Prosthesis
- Durability: Long-lasting
- Precision: Precise Fit For Optimal Performance
- Hygiene: Easy To Clean And Maintain
- Compatibility: Fits Most Dental Implants
- Versatility: Can Be Used For Both Partial And Full Dentures
- Material: Metal
- Highlight: Cylindrical precision attachment denture, Dental Prosthesis precision attachment denture, Telescoper dental prosthesis crown

Product Description

Telescopic Prosthesis

ADS China Dental Lab presents the Telescopic Prosthesis, designed to provide secure fit, natural aesthetics, and long-lasting durability. Perfect for patients looking for reliable dental solutions.

Advantages

Secure Fit and Stability: Offers enhanced retention and minimizes movement for a comfortable fit.

Natural-Looking Aesthetics: Customizable outer crowns blend seamlessly with natural teeth.

Preservation of Natural Teeth: Supports and preserves remaining natural teeth, promoting long-term oral health.

Improved Chewing Function: Secure fit enhances chewing efficiency, supporting a balanced diet.

Long-Lasting Durability: Manufactured with high-quality materials for reliable, long-term use.

Key Features

Telescopic Design: Inner coping and outer crown fit together precisely, offering excellent retention and stability.

Customizable Aesthetics: Outer crown can be tailored to match the color, shape, and appearance of the patient's natural teeth.

Preservation of Natural Teeth: Helps preserve and support remaining natural teeth, preventing further tooth loss.

Enhanced Chewing Function: Improves chewing efficiency with a secure and stable fit.

Durable and Reliable: Crafted using high-quality materials and precision manufacturing techniques for long-lasting use.

Considerations

Oral Hygiene: Maintaining proper oral hygiene and regular dental visits are crucial for the health of natural teeth and supporting structures.

Patient Evaluation: Individual patient evaluation is necessary to determine the suitability of a Telescopic Prosthesis based on oral health and treatment goals.

Trust ADS China Dental Lab

Trust ADS China Dental Lab for your Telescopic Prosthesis needs. Our commitment to quality and precision ensures each prosthesis provides the perfect blend of stability, aesthetics, and functionality.



+8613430498390

info@chinaadsdentallab.com

adschinadentallab.com

Building 4, Hengguangyao Industrial Park, Yonghe Road, Fuyong, Shenzhen, China 518103